



Glove Products

Protegrity® SLK with Neu-Thera® Sterile Latex Powder-Free Surgical Gloves

- Smooth, anti-slip finish for improved dexterity
- Thinner glove (8.3 mils) provides excellent tactile sensitivity during delicate procedures
- First layer of protection in a double-gloving environment when paired with any Cardinal Health powder-free surgical glove
- Inner coating includes
 - Neu-Thera® formulation to protect, restore, moisturize and soothe your skin to aid in skin wellness and hand hygiene
 - Triple-dip nitrile for effortless donning and extra glove strength
- Low-protein formulation
- Wider palm and shorter fingers for a more comfortable, ergonomic fit when single- or double-gloving
- Interlocking cuff prevents glove rolldown



You and your hands deserve only the best! Pamper yourself with our Neu-Thera® treatment and the silky soft protection layer of our Protegrity® SLK surgical gloves.



Length and Thickness

Minimum length measured from the tip of the middle finger to the cuff.

	Length	Cuff Thickness	Palm Thickness	Finger Thickness
ASTM D3577				
Minimum Limit	10.4 in. / 265mm	3.9 mils / 0.1mm	3.9 mils / 0.1mm	3.9 mils / 0.1mm
Average	11.9 in. / 302mm	6.0 mils / 0.15mm	7.5 mils / 0.19mm	7.8 mils / 0.20mm

Physical Properties

Meet or exceed ASTM D3577, "Standard Specification for Rubber Surgical Gloves".

Unaged	ASTM Limit	Cardinal Health Process Average
Tensile Strength	≥24 MPa	25.5 MPa
Tensile Stress at 500%	≤5.5 MPa	2.3 MPa
Elongation (%)	≥750%	924%

Barrier Protection

All gloves must meet a certain Acceptable Quality Level (AQL) as established by the Food and Drug Administration (FDA). This refers to their freedom from holes and their subsequent level of barrier protection. Gloves with a lower actual AQL will have fewer barrier defects. Cardinal Health's internal requirements are significantly more stringent than FDA and ASTM requirements.

	FDA Surgical (2008)	ASTM D3577	Cardinal Health
AQL Limit	2.5 (1.5)	1.5	1.0

Bacteriophage Penetration

This glove passed ASTM F1671, "Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Bloodborne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System." This method has been specifically defined for modeling the viral penetration of Hepatitis (B and C) and Human Immunodeficiency Viruses transmitted in blood and potentially infectious body fluids. It assesses the effectiveness of materials used in protective clothing for protecting the wearer against contact with bloodborne pathogens using a surrogate microbe suspended in simulated body fluid under conditions of continuous contact. The outcome is either "pass" or "fail." Cardinal Health tests a statistically significant sample size of 32 gloves instead of three as called for in the ASTM method.

Proteins and Allergens

These gloves contain 50 micrograms or less of total water-extractable protein per square decimeter. (Caution: Safe use of these gloves by or on latex-sensitive individuals has not been established.) Total protein content is monitored using ASTM D5712.

Ordering Information

Catalog No.	Size	Catalog No.	Size
2D73TPS55	5½	2D73TPS75	7½
2D73TPS60	6	2D73TPS80	8
2D73TPS65	6½	2D73TPS85	8½
2D73TPS70	7	2D73TPS90	9

Gloves are packaged 50 pairs/box, 4 boxes/case.

Protegrity and Neu-Thera are registered trademarks and *Essential to care* is a trademark of Cardinal Health, Inc. or one of its subsidiaries.

© Copyright 2007. Cardinal Health, Inc. or one of its subsidiaries. All rights reserved.
Lit. No. 2GLV0640 (1207/xM/MC3428)

Cardinal Health
Glove Products
1500 Waukegan Road
McGaw Park, IL 60085

cardinalhealth.com/gloves



CardinalHealth
Essential to care™