


Cardinal Health™ PRO/PRO to GO Negative Pressure Wound Therapy

Patient Quick Reference Guide

About Your Cardinal Health™ NPWT PRO / PRO to GO

Your doctor has chosen the Cardinal Health™ NPWT PRO/PRO to GO to remove fluid from your wound by using carefully controlled suction. It is important, however, for you to carefully watch the wound and the PRO/PRO to GO to make sure that it is working properly. Below is some important information and questions that you should ask your healthcare professional.

NOTE: Please refer to the Cardinal Health™ NPWT PRO / PRO to GO PATIENT User Manual for detailed instructions and review the Precautions and Contraindications before use.



Things you need to know about your PRO/PRO to GO

- Do not allow the PRO/PRO to GO to get wet. Clamp the tubing and disconnect from the canister if you take a bath or shower.
- Keep the PRO/PRO to GO plugged in whenever possible to keep the battery fully charged. Always take the A.C. Power Adapter with you when you leave home.
- Keep the PRO/PRO to GO upright to avoid a false Canister Full alert.
- Keep the PRO/PRO to GO turned on at all times unless there is bleeding from the wound or instructed by your healthcare professional.
- Do not change the settings on the PRO/PRO to GO unless you are told to do so by your healthcare professional.

Things to ask your Healthcare Professional

- How to tell if there is a problem with your PRO/PRO to GO or dressing.
- What to do if you have a problem or a leak with your dressing.
- What to do if you notice bleeding from the wound.
- What to do if you must take your dressing off.
- What activities you can do while using the PRO/PRO to GO.
- Who to call if you need help.
- How to take care of your PRO/PRO to GO.

Changing the Canister

1. Clamp the tubing closed.
2. Turn the PRO/PRO to GO off by pressing the OFF **OFF** Button.
3. Remove the blue tapered connector from the top of the canister. Twisting the blue tapered connector makes removing the tubing from the canister easier.
4. To remove the canister, press the Canister Release Button located above the canister, grasp the canister at the bottom and pull downward.
5. Cap and dispose of the used canister properly. Ask your healthcare professional how to properly dispose of a used canister.
6. To install a new canister, hold the new, unused canister at the bottom and slide upwards into the Canister Holder.
7. Align the two short ports and press until it “clicks” into place.
8. Reconnect the tubing to the canister. Gently twist and push the blue tapered connector on the open port of the canister just enough to secure and seal it.
9. Reopen the tubing clamp.
10. Turn the PRO/PRO to GO ON by pressing the ON **ON** Button.

Troubleshooting

If the PRO/PRO to GO begins to alert, refer to the troubleshooting table on the back.

Turning the PRO/PRO to GO ON

Press the ON **ON** Button to start therapy.

Turning the PRO/PRO to GO OFF

To turn PRO/PRO to GO off, press the OFF **OFF** Button. The Alert Display will flash “o-f-f” three times while the PRO/PRO to GO is preparing to shut down.


For questions or information, contact Cardinal Health at 1.866.484.6798

Notice to Users: As with any prescription medical device, failure to follow product instructions or adjusting settings and performing therapy applications without the express direction and/or supervision of your trained clinical caregiver may lead to improper product performance and the potential for serious or fatal injury.

Defibrillation: Remove the wound dressing if defibrillation is required in the area of dressing placement. Failure to remove the dressing may inhibit electrical current transmission and/or patient resuscitation.

Magnetic Resonance Imaging (MRI): The PRO/PRO to GO is not MRI-compatible. Do not take into the MRI area.

Hyperbaric Oxygen Therapy (HBO): NEVER allow a PRO/PRO to GO — whether on or off — inside a hyperbaric chamber. The device must be disconnected from the patient prior to HBO treatment.

What you see or hear	Problem	What to do	More Information
<p>Low Pressure/Leak "1" flashing in Alert Display. Single beep. PRO/PRO to GO is making more noise.</p>	<p>There is an air leak in either the dressing or the tubing connections.</p> <p>Leaks often occur over areas of moist skin, creases or folds in skin, and wrinkles in the drape. They can occur if the drape snags on clothing or bedsheets.</p>	<ul style="list-style-type: none"> • Clamp the tubing. • Look for leaks in the dressing. Gently press around the dressing to check for leaks. If leak is found, use the extra drape or drape strips to seal. • If the alert continues, check the tubing connection to the canister and make sure the tubing is secure by twisting the connector clockwise. • If alert continues, make sure the canister is fully seated and locked in the PRO/PRO to GO. Check for cracks in the canister. If found, replace the canister. • If alert continues, check the tubing for leaks. If a leak is found in the tubing, the entire dressing and tubing set must be replaced. Contact your healthcare professional. • Open the clamp. 	<p>If the leak is properly sealed, the pump becomes quiet and the alert stops. If alert continues, call Cardinal Health at 1.866.484.6798.</p>
<p>Canister Full "2" flashing in Alert Display. Two beeps.</p>	<p>The canister is full.</p>	<ul style="list-style-type: none"> • Make sure clamp is open. • Visually assess the canister to see if full. If the canister is full, change the canister. • If the canister is not full, turn the PRO/PRO to GO off by pressing the OFF  button and then turn the PRO/PRO to GO back on to resume therapy. 	<p>The Canister Full alert begins when the canister is 90 percent full, but the PRO/PRO to GO continues to work until the canister completely fills.</p> <p>If the PRO/PRO to GO is placed on its front, fluid entering the canister will cause a false Canister Full alert and the canister must be changed.</p> <p>If alert continues, call Cardinal Health at 1.866.484.6798 for more assistance.</p>
<p>Low Battery "3" flashing in Alert Display. Three beeps.</p>	<p>The battery is low and has approximately 30 minutes before the battery will be too low to support continued operation of the PRO/PRO to GO.</p>	<p>Plug in the PRO/PRO to GO. A green light shows next to the Plugged In symbol and a yellow light shows next to the Battery Charging symbol to indicate that the battery is charging. The yellow light turns off after battery is fully charged.</p>	<p>Use only the A.C. Power Adapter that came with the PRO/PRO to GO.</p> <p>If alert continues or replacement A.C. Power Adapter is needed, call Cardinal Health at 1.866.484.6798 for assistance.</p>
<p>Service Needed "4" flashing in Alert Display. Four beeps every 10 seconds.</p>	<p>PRO/PRO to Go needs to be serviced.</p>	<p>Call your healthcare professional.</p>	<p>This alert cannot be muted or manually reset by turning the PRO/PRO to GO off and on.</p>
<p>Pressure Setting will not change.</p>	<p>Pressure Setting lock-out is engaged.</p>		<p>Your Healthcare Professional has locked your device to prevent accidental therapy changes.</p> <p>Ask your Healthcare Professional for help.</p>
<p>PRO/PRO to GO is quiet and fluid is not moving in the tubing.</p>	<p>This is NOT a problem.</p>	<p>No action needed.</p>	<p>When the dressing has a good seal, fluid may be removed from the wound and stay in the tubing. The foam is compressed normally and the PRO/PRO to GO is quiet.</p>
<p>An amber light is showing on the front of the PRO/PRO to GO below the pressure numbers.</p>	<p>This is NOT a problem. The PRO/PRO to GO is operating in Intermittent Mode.</p>	<p>No action needed.</p>	<p>Intermittent Mode maintains target pressure for 5 minutes and decreases to -25mmHg for 2 minutes.</p>