



**Presource® Procedure Pack Solutions**

# Anesthesia Kits



**CardinalHealth**



*Comprehensive  
anesthesia solutions*

Safety focus

Procedural efficiencies

Patient comfort



**CardinalHealth**

# Safety, Comfort, Convenience

We offer a new line of standard anesthesia kits that feature quality, latex-free, single-use components. In addition to our selection of standard kits, we can customize a kit to meet your needs using our extensive catalog of components.

By combining all anesthesia supplies in one convenient kit, you are able to achieve efficiencies that enable you to focus on the safety and comfort of the patient.

## **Our anesthesia kits provide:**

- Versatility in the OR through multiple mask and circuit configurations
- Assurance that safe, high-quality components are used in every procedure
- Standardization of anesthesia products and setup
- Reduced on-hand inventory and waste
- Efficiencies in ordering, receiving, case pick and setup
- Improved ease of charge capture

## **Standard or custom options**

Our anesthesia kits come in both standard and custom configurations. If you have specific needs that are not met with our standard kit offering, please contact your **Presource**<sup>®</sup> sales representative to create a custom kit.

# Standard Product Offering

## ANE90ANZ01

### Adult Anesthesia Kit With HME Filter

5/cs

- 1 Gas sampling line, 10' M/M
- 1 Filtered HME
- 1 33" to 87" expandable dual-limb circuit, 3L latex-free bag
- 1 BioTac® electrodes, 5/pk
- 1 Suction tubing, 6mm x 9'
- 1 Yankauer suction handle
- 1 Premium Plus™ mask, regular adult



## ANE90ANZ05

### Adult Anesthesia Kit

5/cs

- 1 Gas sampling line, 10' M/M
- 1 33" to 87" expandable dual-limb circuit, BV filter, 3L latex-free bag
- 1 BioTac® electrodes, 5/pk
- 1 Suction tubing, 6mm x 9'
- 1 Yankauer suction handle
- 1 Premium Soft Plus™ mask, regular adult



## ANE90ANZ02

### Adult Anesthesia Kit With HME and Head Positioner

5/cs

- 1 Gas sampling line, 10' M/M
- 1 Filtered HME
- 1 33" to 87" expandable dual-limb circuit, 3L latex-free bag
- 1 BioTac® electrodes, 5/pk
- 1 Suction tubing, 6mm x 9'
- 1 Yankauer suction handle
- 1 Foam cradle head positioner
- 1 Premium Soft Plus™ mask, regular adult



## ANE90ANZ06

### Adult Anesthesia Kit

5/cs

- 1 Gas sampling line, 10' M/M
- 1 33" to 87" expandable dual-limb circuit, BV filter, 3L latex-free bag
- 1 BioTac® electrodes, 3/pk
- 1 Premium™ mask, regular adult
- 1 Suction tubing, 6mm x 9'
- 1 Yankauer suction handle



## ANE90ANZ03

### Adult Anesthesia Kit With Head Positioner

5/cs

- 1 Gas sampling line, 10' M/M
- 1 33" to 87" expandable dual-limb circuit, BV filter, 3L latex-free bag
- 1 Premium Soft Plus™ mask, regular adult
- 1 BioTac® electrodes, 5/pk
- 1 Suction tubing, 6mm x 9'
- 1 Yankauer suction handle
- 1 Foam cradle head positioner



## ANE90ANZ07

### Adult Anesthesia Kit, Single Limb

5/cs

- 1 Gas sampling line, 10' M/M
- 1 72" single-limb circuit, insp. filter, 3L latex-free bag
- 1 BioTac® electrodes, 5/pk
- 1 Suction tubing, 6mm x 9'
- 1 Yankauer suction handle
- 1 Premium Soft Plus™ mask, regular adult



## PZ90LFANG

### Adult Anesthesia Kit

5/cs

- 2 Exam gloves, powder-free, sterile
- 1 40" dual-limb circuit, BV filter, 3L latex-free bag
- 1 Premium Soft Plus™ mask, regular adult
- 1 Suction tubing, 6mm x 9'
- 1 Yankauer suction handle
- 1 BioTac® electrodes, 3/pk
- 2 3cc luer lock syringe with 21G x 1" needle
- 3 10cc Luer lock syringe with 21G x 1½" needle
- 1 Suction catheter, 18 Fr



## ANE90ANZ04

### Pediatric Anesthesia Kit

5/cs

- 1 Gas sampling line, 10' M/M
- 1 Pediatric 40" to 100" expandable dual-limb circuit, BV filter, 1L latex-free bag
- 1 Medi-Trace® pediatric electrodes, 5/pk
- 1 Premium Plus™ mask, child
- 1 Suction tubing, 6mm x 9'
- 1 Yankauer suction handle



# Suggested Components

ANE90ANZ01  
ANE90ANZ02  
ANE90ANZ03  
ANE90ANZ04  
ANE90ANZ05  
ANE90ANZ06  
ANE90ANZ07

Qty.	Description	Filter	Breathing Bag	Vendor	ANE90ANZ01	ANE90ANZ02	ANE90ANZ03	ANE90ANZ04	ANE90ANZ05	ANE90ANZ06	ANE90ANZ07	
<b>Circuits*</b>												
<b>Adult</b>	_____	Dual, expandable, 33" to 87"	BV <sup>+</sup>	3L	Smiths Medical			1		1	1	
	_____	Dual, expandable, 33" to 87"		3L		1	1					
	_____	Dual, 40"	BV <sup>+</sup>	3L								
	_____	Dual, 60"	BV <sup>+</sup>	3L								
	_____	Single, 40"	Insp. <sup>++</sup>	3L								
	_____	Single, 72"	Insp. <sup>++</sup>	3L								1
<b>Pediatric</b>	_____	Dual, 40"	BV <sup>+</sup>	1L								
	_____	Dual, expandable, 40" to 100"	BV <sup>+</sup>	1L			1					
	_____	Single, 40"	Insp. <sup>++</sup>	2L								
	_____	Single, 60"	Insp. <sup>++</sup>	1L								
<b>Masks</b>												
<b>Premium™</b>	_____	Regular adult			Smiths Medical						1	
	_____	Child										
<b>Premium Plus™</b>	_____	Regular adult with inflation valve			Smiths Medical	1			1		1	
	_____	Child with inflation valve						1				
<b>Premium Soft Plus™</b>	_____	Regular adult with inflation valve and extra cushion			Smiths Medical		1	1				
	_____	Child with inflation valve and extra cushion										
<b>Additional Components</b>												
<b>Filters</b>	_____	Bacterial viral filter			Smiths Medical							
	_____	Filtered heat and moisture exchange (HME)				1	1					
	_____	Bacterial viral filter with HEPA										
<b>Gas Sampling Lines</b>	_____	10" x .050" ID, male/male			Smiths Medical	1	1	1	1	1	1	
	_____	10" x .050" ID, male/female										
<b>Suction Handles and Tubing</b>	_____	Yankauer, bulbous tip, vented			Cardinal Health	1	1	1	1	1	1	
	_____	Yankauer, bulbous tip, nonvented				1	1	1	1	1	1	
	_____	Suction tubing, 6mm x 9', nonsterile										
<b>Electrodes</b>	_____	BioTac® ultra electrodes, 5/pk			Tyco/Kendall	1	1	1		1	1	
	_____	BioTac® ultra electrodes, 3/pk									1	
	_____	Medi-Trace® pediatric electrodes, 5/pk						1				
<b>Head Positioner</b>	_____	Foam positioning cradle, adult			Cardinal Health/Preferred	1	1					
<b>Syringe/Needle Combinations</b>	_____	3cc Luer lock syringe with 21G x 1½" BD Precision Glide™ needle			BD Medical							
	_____	5cc Luer lock syringe with 21G x 1½" BD Precision Glide™ needle										
	_____	10cc Luer lock syringe with 20G x 1½" BD Precision Glide™ needle										

\*All circuits have a GSP at the elbow. All dual-limb circuits have a straight wye. +BV = Bacterial Viral, ++Insp. = Inspiratory

## General Customer Information

Catalog Number	Sell Price (each)	
Customer Name	Date	Monthly Usage
Contact Name	Contact Title	Contact Phone
Special Instructions		

## Presource® Representative Use Only

SOR \_\_\_\_\_ Sample Required  No  Yes, Date Requested \_\_\_\_\_

Presource® Representative \_\_\_\_\_ Customer Number \_\_\_\_\_ Contact Phone \_\_\_\_\_ Date Quote Requested \_\_\_\_\_


Please fax this page to 614.652.7268.

Fold and tear on dotted line before faxing.

## Presource® Procedure Pack Solutions

**Presource** is a registered trademark and *Essential to care* is a trademark of Cardinal Health, Inc. or one of its subsidiaries. Premium, Premium Plus and Premium Soft Plus are trademarks of Portex Corporation. BioTac and Medi-Trace are registered trademarks of Tyco Healthcare Group LP. BD PrecisionGlide is a trademark of Becton, Dickinson and Company.

© Copyright 2008. Cardinal Health, Inc. or one of its subsidiaries. All rights reserved.  
Lit. No. 2PRE0791-01 (0708/xM/MC3663)

 Paper contains a minimum of 10% post-consumer fiber.



Cardinal Health  
**Presource®** Products  
7000 Cardinal Place  
Dublin, OH 43017

[cardinalhealth.com/presource](http://cardinalhealth.com/presource)