



Presource® Procedure Pack Solutions

Dressing Change Kits




CardinalHealth

Stay ahead of infection by promoting a safe, consistent standard of care

Reducing Hospital-Acquired Infections (HAIs) is imperative to reducing healthcare costs and to improving quality of care and quality of life for patients.

- The Centers for Disease Control (CDC) states that 250 people die every day from HAIs.
- HAIs erode profits, adding \$4.5 to \$5.7 billion dollars in healthcare costs.

Annually, bloodstream infections alone account for more than 248,000 infections and 30,665 deaths.



The staggering number of infections has led to the development of new guidelines for standards of care and an increase in government intervention.

As knowledge, technology and healthcare settings change, infection control and prevention measures should also change.

Cardinal Health offers these products which can help prevent infection

ChloroPrep® Solution

The 2% Chlorhexidine Gluconate (CHG) and 70% Isopropyl Alcohol (IPA) formulation acts fast on a broad spectrum of bacteria, including Methicillin-Resistant Staphylococcus aureus (MRSA), and keeps fighting bacteria for at least 48 hours.¹

Compared to Povidone Iodine Solution, CHG has shown a 50% reduction in catheter-related bloodstream infections.

The CDC Guidelines state that disinfecting skin with 2% Chlorhexidine Gluconate is preferred at time of dressing change.



BIOPATCH® Protective Disk With Chlorhexidine Gluconate (CHG)

BIOPATCH® delivers powerful skin antisepsis for up to a 7-day time period in between dressing changes, creating a zone of inhibition for the surrounding area of the site.



Research shows that using BIOPATCH® on patients with central venous and arterial catheters reduced the incidence of local infections by 44% and bloodstream infections by 60% as compared to standard care.

Suggested Components

03-3000A
03-3200A
NI03-0130B
03-4500A
24012-060
26012-633
26012-639
3T3008
NI03-0515B
NI03-0544A
03-2000
03-2100

Qty.	Description	Vendor	03-3000A	03-3200A	NI03-0130B	03-4500A	24012-060	26012-633	26012-639	3T3008	NI03-0515B	NI03-0544A	03-2000	03-2100
Drape	Drape, 13" x 19"	Cardinal Health/ Preferred			1	1		1						1
	Drape, 17" x 19"	Cardinal Health/ Preferred								1				
Dressing	Tegaderm® dressing, 3½" x 4½"	3M			1									1
	Tegaderm® dressing, 4" x 4¾"	3M	1			1		1	1	1	1	1	1	
	IV split dressing, 2" x 2"	Cardinal Health/ Preferred					1							
	IV stabilizing dressing, 4" x 5½"	Smith & Nephew		1			1							
	Dressing, 2" x 3"	Cardinal Health/ Preferred	1	1		1			1	1				
	Tegaderm® dressing, IV film	3M												
Gauze	BIOPATCH® protective disc, 1"	Ethicon											1	1
	Gauze, 2" x 2", 8-ply	Cardinal Health/ Preferred	1	1		1			1	1	3	2	3	
	Gauze, 3" x 3", 12-ply	Cardinal Health/ Preferred												
	Gauze, 4" x 4", 8-ply	Cardinal Health/ Preferred	1	1	2	1	2		1	1				2
	Gauze sponge, 2" x 2", 4-ply	Kendall												
Glove	Sponge dressing, 2" x 2", 6-ply	Kendall	1	1		1		1	1					
	Pair of gloves, vinyl, powder-free, cuffed	Cardinal Health/ Preferred										2		
	Pair of gloves, vinyl, powdered, cuffed	Cardinal Health/ Preferred						2	2					
	Pair of gloves, powder-free, cuffed, walled	Cardinal Health/ Preferred	1	1	1	1	1			1	1		1	
Instrument	Pair of gloves, powder-free, walled	Cardinal Health/ Preferred												1
	Scissors, S/B, 4½" wire	Cardinal Health/ Preferred					1							
Label	Label, change date and initials, 1" x 2"	Cardinal Health/ Preferred	1		1		1			1				1
Mask	Mask	Cardinal Health/ Preferred	1	1	1	1			1		1	1	1	1
Misc.	Refuse bag, 6" x 12" x 18"	Cardinal Health/ Preferred	1											
	Refuse bag with tie, 12" x 18" x 6"	Cardinal Health/ Preferred		1			1							
Prep	PVP ointment	Cardinal Health/ Preferred		1		1	1	1	1	1				
	Triple antibiotic ointment	Cardinal Health/ Preferred	1											
	Acetone/alcohol prep pad	Cardinal Health/ Preferred					3							
	Alcohol prep pad	Cardinal Health/ Preferred	3	3										
	Skin protectant prep pad	Cardinal Health/ Preferred	1	1			1							
	ChloroPrep® one-step prep stick, 3mL	Cardinal Health				1					1	1	1	1
	PVP swab stick triples	Cardinal Health	1	1			1	1	1	1				
	Acetone/alcohol swab stick	Cardinal Health					1							
	Alcohol swab stick triples	Cardinal Health	1	1			1	1	1	1				
	Skin protectant swab stick	Cardinal Health				1		1						
Tape	Plastic applicator, 3"	Cardinal Health	1	1			1							
	Plastic swab stick, 6"	Cardinal Health					1							
	Tape measure	Cardinal Health											1	
	Transpore® tape, 1" x 18"	3M				1		1	1	1	1	1	1	1
Towel	Transpore® tape, 1" x 36"	3M	1	1		1								
	Plasti-pore® tape, 1" x 18"	Cardinal Health												
Towel	Towel, 17" x 19", paper	Cardinal Health										1		

General Customer Information

Catalog Number _____ Sell Price (each) _____

Customer Name _____ Date _____ Monthly Usage _____

Contact Name _____ Contact Title _____ Contact Phone _____

Special Instructions _____

Presource® Representative Use Only

SOR _____ Sample Required No Yes, Date Requested _____

Presource® Representative _____ Customer Number _____ Contact Phone _____ Date Quote Requested _____

Please fax this page to 614.652.7268.

Fold and tear on dotted line before faxing.

Fold and tear on dotted line before faxing.

Standard Product Offering

26012-639

Dressing Change Kit With Tegaderm® Dressing

30/cs

- 1 Tegaderm® dressing, 4" x 4¾"
- 1 Dressing, 2" x 3"
- 1 Gauze, 2" x 2", 8-ply
- 1 Gauze, 4" x 4", 8-ply
- 1 Sponge dressing, 2" x 2", 6-ply
- 2 Pairs of gloves, vinyl, powdered, cuffed
- 1 Mask
- 1 PVP ointment
- 1 Pack of PVP swab stick triples
- 1 Pack of alcohol swab stick triples
- 1 Transpore® tape, 1" x 18"
- 1 CSR wrap



NI03-0544

Central Line Dressing Change Kit With ChloraPrep® solution and Tegaderm® Dressing

20/cs

- 1 Tegaderm® dressing, 4" x 4¾"
- 2 Gauze, 2" x 2", 8-ply
- 2 Pairs of gloves, vinyl, powder-free, cuffed
- 1 Mask
- 1 ChloraPrep® solution, 3mL
- 1 Infant tape measure, 36"
- 1 Transpore® tape, 1" x 18"
- 1 Towel, 17" x 19"
- 1 CSR wrap



3T3008

Central Line Dressing Change Kit With Tegaderm® Dressing

20/cs

- 1 Drape, 17" x 19"
- 1 Tegaderm® dressing, 4" x 4¾"
- 1 Dressing, 2" x 3"
- 1 Gauze, 2" x 2", 8-ply
- 1 Gauze, 4" x 4", 8-ply
- 1 SORB-IT® sponge, 2" x 2"
- 1 Pair of gloves, powder-free, cuffed, walled
- 1 Label, change date and initials, 1" x 2"
- 1 PVP ointment
- 1 Pack of PVP swab stick triples
- 1 Pack of alcohol swab stick triples
- 1 Transpore® tape, 1" x 18"

03-2000

Central Line Dressing Change Kit With BIOPATCH® Protective Disc

30/cs

- 1 Tegaderm® dressing, 4" x 4¾"
- 1 BIOPATCH® protective disc, 1"
- 3 Gauze, 2" x 2", 8-ply
- 1 Pair of gloves, powder-free, cuffed, walled
- 1 Mask
- 1 ChloraPrep® solution, 3mL
- 1 Transpore® tape, 1" x 18"
- 1 CSR wrap

NI03-0515B

Central Line Dressing Change Kit With ChloraPrep® Solution and Tegaderm® Dressing

30/cs

- 1 Tegaderm® dressing, 4" x 4¾"
- 3 Gauze, 2" x 2", 8-ply
- 1 Pair of gloves, powder-free, cuffed, walled
- 1 Mask
- 1 ChloraPrep® solution, 3mL
- 1 Transpore® tape, 1" x 18"
- 1 CSR wrap



03-2100

Central Line Dressing Change Kit With BIOPATCH® Protective Disc and Drape

30/cs

- 1 Drape, 13" x 19"
- 1 Tegaderm® dressing, 3½" x 4½"
- 1 BIOPATCH® protective disc, 1"
- 2 Gauze, 4" x 4", 8-ply
- 1 Pair of gloves, powder-free, walled
- 1 Label, change date and initials, 1" x 2"
- 1 Mask
- 1 ChloraPrep® solution, 3mL
- 1 Transpore® tape, 1" x 18"
- 1 CSR wrap

Standard Product Offering *(continued)*

03-3000A

Dressing Change Kit With Tegaderm® Dressing

30/cs

- 1 Tegaderm® dressing, 4" x 4¾"
- 1 Dressing, 2" x 3"
- 1 Gauze, 2" x 3", 8-ply
- 1 Gauze, 4" x 4", 8-ply
- 1 Sponge dressing, 2" x 2", 6-ply
- 1 Pair of gloves, powder-free, cuffed, walled
- 1 Label, change date and initials, 1" x 2"
- 1 Mask
- 1 Refuse bag, 6" x 12" x 18"
- 1 Triple antibiotic ointment
- 3 Alcohol prep pads
- 1 Skin protectant prep pad
- 1 Pack of PVP swab stick triples
- 1 Pack of alcohol swab stick triples
- 1 Plastic applicator, 3"
- 1 Transpore® tape, 1" x 36"
- 1 CSR wrap

03-3200A

Dressing Change Kit With Dressing

30/cs

- 1 IV stabilizing dressing, 4" x 5½"
- 1 Dressing, 2" x 3"
- 1 Gauze, 2" x 2", 8-ply
- 1 Gauze, 4" x 4", 8-ply
- 1 Sponge dressing, 2" x 2", 6-ply
- 1 Pair of gloves, powder-free, cuffed, walled
- 1 Mask
- 1 Refuse bag, white, 12" x 18" x 6"
- 1 PVP ointment
- 3 Alcohol prep pads
- 1 Skin protectant prep pad
- 1 Pack of PVP swab stick triples
- 1 Pack of alcohol swab stick triples
- 1 Plastic applicator, 3"
- 1 Transpore® tape, 1" x 36"
- 1 CSR wrap

NI03-0130B

Central Line Dressing Change Kit With Chloraprep® Solution and Tegaderm® Dressing

30/cs

- 1 Drape, 13" x 19"
- 1 Tegaderm® dressing, 3½" x 4½"
- 2 Gauze, 4" x 4", 8-ply
- 1 Pair of gloves, powder-free, cuffed, walled
- 1 Mask
- 1 Label, change date and initials, 1" x 2"
- 1 Chloraprep® solution, 3mL
- 1 Transpore® tape, 1" x 18"
- 1 CSR wap

03-4500A

Central Line Dressing Change Kit With Tegaderm® Dressing

30/cs

- 1 Drape, 13" x 19"
- 1 Tegaderm® dressing, 4" x 4¾"
- 1 Dressing, 2" x 3"
- 1 Gauze, 2" x 2", 8-ply
- 1 Gauze, 4" x 4", 8-ply
- 1 Sponge dressing, 2" x 2", 6-ply
- 1 Pair of gloves, powder-free, cuffed, walled
- 1 Mask
- 1 PVP ointment
- 1 Pack of PVP swab stick triples
- 1 Pack of alcohol swab stick triples
- 1 Skin protectant swab stick
- 1 Transpore® tape, 1" x 36"
- 1 CSR wrap

24012-060

Dressing Change Kit With Dressing

20/cs

- 1 IV split dressing, 2" x 2"
- 1 IV stabilizing dressing, 4" x 5½"
- 2 Gauze, 4" x 4", 8-ply
- 1 Pair of gloves, powder-free, cuffed, walled
- 1 Scissors, S/B, 4½", wire
- 1 Label, change date and initials, 1" x 2"
- 1 Refuse bag, white, 12" x 18" x 6"
- 1 PVP ointment
- 3 Acetone/alcohol prep pads
- 1 Skin protectant prep pad
- 1 Pack of PVP swab stick triples
- 1 Pack of acetone/alcohol swab stick triples
- 1 Plastic applicator, 3"
- 1 CSR wrap

26012-633

Dressing Change Kit With Tegaderm® Dressing

50/cs

- 1 Drape, 13" x 19"
- 1 Tegaderm® dressing, 4" x 4¾"
- 1 Sponge dressing, 2" x 2", 6-ply
- 2 Pairs of gloves, powdered, cuffed
- 1 PVP ointment
- 1 Pack of PVP swab stick triples
- 1 Pack of alcohol swab stick triples
- 1 Skin protectant swab stick
- 1 Transpore® tape, 1" x 18"



Presource® Procedure Pack Solutions

¹Garcia R., Mulberry G, Brady A, Hibbard JS. Comparison of ChloroPrep® and BETADINE® as preoperative skin preparation antiseptics. Poster presented at: 40th Annual Meeting of the Infectious Disease Society of America: October 25, 2002; Chicago, IL.

Presource and ChloroPrep are registered trademarks and *Essential to care* is a trademark of Cardinal Health, Inc. or one of its subsidiaries.

BETADINE is a registered trademark of The Perdue Frederick Company.

BIOPATCH is a registered trademark of Ethicon, Inc.

Plasti-pore is a registered trademark of Smith & Nephew, Inc.


SORB-IT is a registered trademark of The Kendall Company LP.

Tegaderm and Transpore are registered trademarks of 3M Company.

© Copyright 2008. Cardinal Health, Inc. or one of its subsidiaries.

All rights reserved.

Lit. No. 2PRE0724 (0608/1M/MC3532)

 Paper contains a minimum of 10% post-consumer fiber.



Cardinal Health
Presource® Products
7000 Cardinal Place
Dublin, OH 43017

cardinalhealth.com/presource