

Cardinal Health™ NPWT Disposables and Accessories

We offer a variety of disposables and accessories for all of your Cardinal Health™ NPWT devices.

Wound dressings

- Black foam kit
Includes:
 - Black foam — open cell, reticulated, hydrophobic, high tensile strength
(Available in small, medium, large and extra-large sizes)
 - Polyurethane drape — easy to apply, semi-permeable for breathability
 - SpeedConnect™ tubing — low-profile flange, pliable, designed to resist kinking
- White foam — hydrophilic, high density, small pore size, flexible when wet or dry, no need to wet the foam prior to placement
(Available in standard and large sizes)

Ordering information

| Item | Size | Cat. No. | Qty. |
|--------------------------|-------------------------|--------------------------|----------|
| Black foam dressing kit | Small, 10 x 8 x 3cm | 47-1702 | 10 ea/cs |
| | | 47-1702/15 | 15 ea/cs |
| | Medium, 20 x 12.5 x 3cm | 47-1701 | 10 ea/cs |
| | | 47-1701/15 | 15 ea/cs |
| | Large, 25 x 15 x 3cm | 47-1700 | 10 ea/cs |
| | | 47-1700/15 | 15 ea/cs |
| XLarge, 58.5 x 33 x 3cm | 47-1703 | 10 ea/cs | |
| | White foam | Small, 9.5 x 9.5 x 0.5cm | 47-1751 |
| Large, 35 x 25 x 0.635cm | | 47-1755 | 5 ea/cs |



Canisters

- Robust, impact resistant design
- Filters to help prevent leaks and control odor
- Available with and without gel packs to solidify exudate
- Easy to install with one snap
- Integrated volume graduation



Ordering information

| Item | Description | Cat. No. | Qty. |
|---------------------------------|-------------|----------|----------|
| 300cc canister | With gel | 47-4000 | 10 ea/cs |
| | Without gel | 47-4100 | 10 ea/cs |
| 500cc canister | With gel | 47-4500 | 10 ea/cs |
| 1200cc canister (Svedman® only) | With gel | 47-3000 | 10 ea/cs |

Accessories

- Bridging kit and Y-connector for re-positioning of the SpeedConnect™ for patient comfort
- Irrigation tubing with SpeedConnect™ tubing and Irrigation delivery kit for delivery of Simultaneous Irrigation™ Technology
- SensiSkin™ drape for patients with sensitive skin
- IV pole adapter and bed brackets for secure and easy placement of devices



Ordering information

| Item | Description | Cat. No. | Qty. |
|---|------------------|----------|----------|
| Polyurethane drape | | 47-7000 | 10 ea/cs |
| SensiSkin™ drape | | 47-7100 | 10 ea/cs |
| Bridging kit | | 47-1704 | 5 ea/cs |
| Y-connector | | 47-2500 | 10 ea/cs |
| SpeedConnect™ tubing set | | 47-2000 | 10 ea/cs |
| Irrigation tubing with SpeedConnect™ tubing set | | 47-6000 | 5 ea/cs |
| Irrigation delivery kit | | 47-6500 | 5 ea/cs |
| PRO device accessories | Carrying case | 47-9600 | 1 ea/cs |
| | IV pole adapter | 47-5600 | 1 ea/cs |
| | AC power adapter | 47-9100 | 1 ea/cs |
| SVED® device accessories | Carrying case | 47-9500 | 1 ea/cs |
| | IV pole adapter | 47-5500 | 1 ea/cs |
| | AC power adapter | 47-9000 | 1 ea/cs |

For more information or to order, contact your Cardinal Health sales representative today, call **866.484.6798** or visit cardinalhealth.com/npwt

