

Cleanroom Gloves

reference guide



Advancing performance by design.

Whether latex or synthetic, sterile or non-sterile, hand-specific or ambidextrous, Cardinal Health offers a portfolio of gloves that will help meet a wide variety of needs. Cardinal Health™ Cleanroom Gloves are clean-processed and designed to meet the demands of pharmaceutical, medical device and biotech manufacturers in controlled environments, as well as professionals in lab research environments.

Sterile

Glove	Cat. no.	ISO	Material	Description
CR100 Sterile Nitrile	2Y1860T – 2Y1867T	4	Nitrile	Double chlorinated to ensure low particulate and extractable counts. Hand specific. Beaded cuff.
CR100 Sterile Latex	2Y1624 – 2Y1630, 2Y1640	4	Latex	Double chlorinated to ensure low particulate and extractable counts. Hand specific
Duraprene™ CP	2Y72PN55 – 2Y72PN90	5	Neoprene	Offers the fit, feel and performance of a latex glove without the latex proteins, and utilizes a nitrile coating help to improve donning and durability. Beaded cuff.
CP100 BT	2Y5000 – 2Y5007	5	Latex	Features a textured, slip-resistant finish. Hand specific. Beaded cuff.
Esteem® CP	2Y72PL55 – 2Y72PL90	6	Polyisoprene	Fit, feel and performance that is virtually indistinguishable from latex. Hand specific. Beaded cuff.
Protegrity® CP	2Y72N0 – 2Y72N7	6	Latex	Unique bonding technology that integrates the comfort of latex with the strength of nitrile coating for grip. Hand specific. Beaded cuff.

Non-Sterile

Glove	Cat. no.	ISO	Material	Description
CR100 Non-Sterile Latex	2Y1616 – 2Y1622, 2Y1639	4	Latex	Double chlorinated to ensure low particulate and extractable counts. Hand specific
Multi-Flex® CC100	2Y1600 – 2Y1607	5	Latex	More flexible fit. Hand specific
CP100 Nitrile Ambi	2Y1849 – 2Y1853	6	Nitrile	Features a textured, split-resistant finish. Ambidextrous. Beaded cuff.
CP100 BT Ambi	2Y1730 – 2Y1733	6	Latex	Features a textured, slip-resistant finish. Ambidextrous. Beaded cuff.
Esteem® Stretchy Nitrile	2Y1870T – 2Y1874T	7	Nitrile	Features a textured, slip-resistant finish. Ambidextrous. Beaded cuff.