


Quick. Accurate. Reliable.

Cardinal Health Clinical Digital Thermometers

Healthcare professionals and their patients need a wide range of temperature monitoring products to ensure quality and accuracy of treatment. From single to dual-scale, our digital thermometers measure temperature rapidly and accurately in either Fahrenheit or Celsius.



Ultra Forehead Temperature Indicator



Rectal thermometer



Oral thermometer



Digital thermometer probe covers

Digital Thermometers

- Easy-to-read LCD display
- Last temperature taken memory capability
- Audible alarm signals, one beep at start-up, 10 beeps max at peak temperature
- Auto turn-off after five (5) minutes of inactivity
- Available in Fahrenheit, Celsius, or switchable Dual Scale
- Standard 60 second read or Quick Read option for faster measurement
- Oral or rectal* thermometers
- Probe covers available separately
- Individually boxed or economical polybags available

*Rectal thermometers come with a red cap for easy identification. Probe covers for rectal thermometers are pre-lubricated.

Ultra Forehead Temperature Indicator

- Encapsulated liquid crystal technology, displays accurate temperature changes in seconds
- Non-invasive and inexpensive alternative to electronic probes
- Indicates temperatures in °C and °F scales with accuracy to $\pm 1^{\circ}$ Fahrenheit.
- Ideal for use with laryngeal masks and mac or regional anesthesia when intubation is not used
- Discrete temperature events
 - A green color indicates the temperature attained
- Not made with natural rubber latex.

For more information or to order, contact your Cardinal Health sales representative, call 877.227.3462 or visit cardinalhealth.com

Cat. no.	O: Oral R: Rectal	F: Fahrenheit C: Celsius DS: Dual Scale	Probe covers /ea	Box or Polybag	Qty.
----------	----------------------	---	------------------	----------------	------

Standard Thermometers (~60 second oral read)

16811ECO	O	F	5	Polybag	25/bx, 500/cs
16811-100	O	F	5	Box	12/bx, 144/cs
16811ECOR	R	F	5	Polybag	25/bx, 500/cs
16811-100R	R	F	5	Box	12/bx, 72/cs
16811-C	O	C	5	Box	12/bx, 144/cs
16811-EDS	O	DS	5	Polybag	25/bx, 500/cs
16811-DS	O	DS	5	Box	12/bx, 144/cs
16811-DSR	R	DS	5	Box	12/bx, 72/cs
16811-500	O	F	0	Polybag	25/bx, 500/cs
16811-500R	R	F	0	Polybag	25/bx, 500/cs


Quick read (~10 second oral read)

16811-EQRO	O	DS	0	Polybag	25/bx, 500/cs
16811-EQRR	R	DS	0	Polybag	25/bx, 500/cs

Probe covers

16811-PRBO	O	n/a	n/a	Box	100/bx, 2000/cs
16811-PRBR	R	n/a	n/a	Box	100/bx, 2000/cs

Cat. no.	Description	Qty.
14677-493	Ultra forehead temperature indicator	100/bx, 2 bx/cs



cardinalhealth.com

Cardinal Health
Medical Products
Waukegan, Illinois