





















## Quick reference









# Cardinal Health™ Exam Gloves



## High use portfolio










A portfolio of non-sterile exam gloves designed to serve a broad range of everyday clinical needs. Highly versatile exam gloves in a variety of styles allow clinicians to select the right product that best meets their protection and performance needs. Select products are also available in premium configurations including features such as dispensing film, increased elongation, low dermatitis potential or moisturizing coating.

Cat. No.										
Fingertip thickness (+/-0.2mil)	Product name	Size	Standard configuration	Premium configuration	Premium feature	Color	Texture/ cuff length	Chemo drug tested	Chemical tested	Packaging
Nitrile										
3.5mil	FLEXAL™ Touch	XS	88RT01XS	88RTF01XS	• Dispensing film	Blue				<b>XS – L:</b> 250 ea/bx; 10 bx/cs  <b>XL – XXL:</b> 230 ea/bx; 10 bx/cs
		S	88RT02S	88RTF02S						
		M	88RT03M	88RTF03M						
		L	88RT04L	88RTF04L						
		XL	88RT05XL	88RTF05XL						
	XXL	88RT06XXL	88RTF06XXL							
3.7mil	FLEXAL™ Comfort	XS		88FS01XS	• Dispensing film • >650% elongation • Low dermatitis	Light blue				<b>XS – L:</b> 250 ea/bx; 10 bx/cs  <b>XL – XXL:</b> 230 ea/bx; 10 bx/cs
		S		88FS02S						
		M		88FS03M						
		L		88FS04L						
		XL		88FS05XL						
		XXL		88FS06XXL						
	FLEXAL™	XS	88TN01XS	88TNF01XS	• Dispensing film	Cornflower blue				<b>XS – XXL:</b> 200 ea/bx; 10 bx/cs
		S	88TN02S	88TNF02S						
		M	88TN03M	88TNF03M						
		L	88TN04L	88TNF04L						
XL		88TN05XL	88TNF05XL							
	XXL	88TN06XXL	88TNF06XXL							
4.5mil	ESTEEM™ Comfort	XS		88ES01XS	• Dispensing film • >650% elongation • Low dermatitis	Light blue				<b>XS – L:</b> 200 ea/bx; 10 bx/cs  <b>XL:</b> 180 ea/bx; 10 bx/cs
		S		88ES02S						
		M		88ES03M						
		L		88ES04L						
		XL		88ES05XL						
	ESTEEM™ Stretch	XS	8854NXSB			Teal blue				<b>XS – L:</b> 150 ea/bx; 10 bx/cs  <b>XL – XXL:</b> 130 ea/bx; 10 bx/cs
		S	8855NSB							
		M	8856NMB							
		L	8857NLB							
		XL	8858NXLB							
	XXL	8859NXXB								
4.7mil	ESTEEM™ Moisturize with Neu-Thera™ Coating	XS		N88RX01T	• Moisturizing coating	Blue				<b>XS – L:</b> 100 ea/bx; 10 bx/cs  <b>XL – XXL:</b> 90 ea/bx; 10 bx/cs
		S		N88RX02T						
		M		N88RX03T						
		L		N88RX04T						
		XL		N88RX05T						
		XXL		N88RX06T						
	ESTEEM™	XS	8895NB			Blue				<b>XS – L:</b> 150 ea/bx; 10 bx/cs  <b>XL – XXL:</b> 130 ea/bx; 10 bx/cs
		S	8896NB							
		M	8897NB							
		L	8898NB							
		XL	8899NB							
	XXL	8899NXXB								

<b>Key</b> USP800 compliant	Thiotepa minimum breakthrough detection time in minutes	Carmustine minimum breakthrough detection time in minutes	Smooth	Fingertip texture	Full finger texture	Finger and palm textured	Fully textured
							

## Specialty portfolio

A portfolio of non-sterile and sterile exam gloves designed to meet specific clinical needs. Exam gloves with features designed to protect clinicians in highly specific applications, such as chemotherapy drug compounding, sterile, high risk or fine tactile sensitivity procedures, and surgical instrument decontamination.

Fingertip thickness (+/-0.2mil)	Specialty feature and/or application	Product name	Size	Cat. No.	Color	Texture/ cuff length	Sterile	Hand- specific	Chemo drug tested	Chemical tested	Packaging
Nitrile											
3.3mil	High tactile sensitivity/ phlebotomy	FLEXAL™ Feel	XS	88TT20XS	Blue					✓	XS – L: 300 ea/bx; 10 bx/cs  XL – XXL: 250 ea/bx; 10 bx/cs
			S	88TT21S							
			M	88TT22M							
			L	88TT23L							
6.2mil	STERILE Sterile procedures	Sterile (singles)	XL	88TT24XL	Blue		✓		✓	✓	S – XL: 50 ea/bx; 5 bx/cs
			S	N8820							
			M	N8821							
			L	N8822							
		Sterile (pairs)	XL	N8823							S – XL: 40 pr/bx; 5 bx/cs
			S	N8830							
M	N8831										
L	N8832										
XL	N8833										
6.9mil	High risk/ chemotherapy drug handling	Extended Cuff	XS	N8850XPB	Blue				✓	✓	XS – XL: 100 ea/bx; 10 bx/cs
			S	N8851XPB							
			M	N8852XPB							
			L	N8853XPB							
8.0mil	Chemotherapy drug handling	CHEMOPLUS™	XL	N8854XPB	Blue				✓	✓	S – XL: 100 ea/bx; 4 bx/cs
			S	CT5071G-							
			M	CT5072G-							
			L	CT5073G-							
	STERILE Sterile chemotherapy drug handling	CHEMOPLUS™ Sterile	XL	CT5074G-	Blue		✓		✓	S – XL: 50 pairs/bx; 4 bx/cs	
			S	CTS731S--							
M	CTS731M--										
L	CTS731L--										
XL	CTS731XL--										
11.0mil	Decontamination	Decontamination	S	88NDS	Blue				✓	✓	S – XL: 50 ea/bx; 10 bx/cs
			M	88NDM							
			L	88NDL							
			XL	88NDXL							
Neoprene											
9.0mil	Chemotherapy drug handling	CHEMOPLUS™	S	CT21911--	Brown			✓	✓		S – XL: 50R and 50L ea/bx; 3 bx/cs
			M	CT21921--							
			L	CT21931-							
			XL	CT21941--							
Latex											
14.1mil	High risk/ decontamination	XP Plus	S	L88HRS	Blue				✓	✓	S – XL: 50 ea/bx; 10 bx/cs
			M	L88HRM							
			L	L88HRL							
			XL	L88HRXL							
18.0mil	Decontamination	CHEMOPLUS™	S	CT0191-1-	Blue			✓	✓		S – XL: 50R and 50L ea/bx; 10 bx/cs
			M	CT0192-1-							
			L	CT0193-1-							
			XL	CT0194-1-							

† Meets the ASTM D 6978 standard

This Cardinal Health product helps healthcare facilities comply with the USP <800> standard for the safe handling of hazardous drugs.

© 2019 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO, ESSENTIAL TO CARE, ESTEEM, FLEXAL, NEU-THERA and CHEMOPLUS are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. Lit. No. 2CHL18-830488 (9/2018)