

Cardinal Health™ Fecal Specimen Transport



- Built in collection spoon for easy sampling
- Easy to read specimen fill line
- Color coded cap and label for ease of use
- Several options to meet your laboratory needs based on your sampling criteria

| Cat. No. | Description | Intended use | Packaging |
|------------------|--|--|-------------------|
| CHB900612 | Culture and Sensitivity stool transport vial | Non-nutritive medium with indicator for enteric pathogen transport | 20 ea/pk, 6 pk/cs |
| CHB900812 | Enteric stool transport vial | Non-nutritive medium with indicator for enteric pathogen transport | 20 ea/pk, 6 pk/cs |
| CHB900412 | 10% Formalin stool transport vial | For ova and parasite stool specimens | 20 ea/pk, 6 pk/cs |
| CHB901012 | Zn-PVA stool transport vial | For ova and parasite stool specimens | 20 ea/pk, 6 pk/cs |
| CHB900212 | SAF stool transport vial | For ova and parasite stool specimens | 20 ea/pk, 6 pk/cs |
| CHB900112 | LV-PVA stool transport vial | For ova and parasite stool specimens | 20 ea/pk, 6 pk/cs |
| CHB901612 | Single Vial stool transport vial | For ova and parasite stool specimens | 20 ea/pk, 6 pk/cs |
| CHB901312 | EcoFix® stool transport vial | For ova and parasite stool specimens | 20 ea/pk, 6 pk/cs |
| CHB900312 | Clean stool transport vial | Empty vial for stool transport | 20 ea/pk, 6 pk/cs |



We are your lab partner.
Expertise. Choice. Efficiency. Savings.

For more information or to order,
contact your Cardinal Health sales representative,
call 888.444.5440 or visit cardinalhealth.com/medical.



EcoFix is a registered trademark of Meridian Biosciences, Inc.

© 2018 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO and ESSENTIAL TO CARE are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. Lit. No. 2LAB18-789758 (4/2018)

cardinalhealth.com