

Specifications

CAHFT30 Bitter (Denatonium Benzoate) Fit Test Kit

Kit:

- Test hood
- Sensitivity nebulizer
- Test nebulizer with atomizer
- Wire cleaner
- Box of sensitivity solution, 6 ampules
- Box of test solution, 6 ampules

Hood:

- Poly-coated
- Clear viewing window
- Surged seams
- Reusable

Bitter (denatonium benzoate):

- Glass ampule
- With pre-drilled plastic cap for easy pouring
- Plastic protection from broken glass
- One-time use
- Averages 3-5 fit tests per tube
- Contains 2.5cc of solution

Nebulizer:

- Buna-N (nitrile rubber) air bulb for easy operation
- Converts 5cc (mL) of liquid into aerosol 0.3-5.0 microns particle size
- Complete with extra jet and replaceable plastic caps
- Break-resistant polycarbonate plastic (Lexan®) that easily disassembles for washing
- Dimensions: 2 in. x 2 in. x 7 in. (5.0 cm x 5.0 cm x 17.7 cm)



This kit offers qualitative OSHA compliant fit testing for disposable and reusable dust and mist respirators. Denatonium benzoate creates an unmistakable bitter taste that indicates breakthrough in the respirator fit. New disposable glass ampules make pouring solution into nebulizers easier and reduces the chance of contamination.

Products available:

CAHFT30	Bitter (Denatonium Benzoate) Fit Test Kit
CAHFT31	Bitter (Denatonium Benzoate) Sensitivity Solution (6 EA/BX)
CAHFT32	Bitter (Denatonium Benzoate) Test Solution (6 EA/BX)
CAHFT10	Nebulizer (1 EA)

For more information or to order, contact your Cardinal Health sales representative or visit [cardinalhealth.com/facialprotection](https://www.cardinalhealth.com/facialprotection).