

Prepare wrinkle-free and clear paraffin embedded tissue sections every time

Cardinal Health Tissue Flotation Baths



Cardinal Health Deluxe Lighted Tissue Flotation Bath

Material code: M7654-5A

- Eliminate wrinkles and distortion during preparation of paraffin embedded tissue sections
- L.E.D. illuminated black background with transparent PYREX® dish; allows specimen to be easily seen
- Digital microprocessor control with advanced safety over temperature protection
- Slide Dryer Block to assist in drying of slides
- Tissue Orienter Block to smooth wrinkled sections

Supplied with:

- Removable PYREX® glass dish (part number: 904-0013)
- Certified to UL61010-1 and CSA C22.2



Cardinal Health Lighted Tissue Flotation Bath

Material code: M7654-1A

- Eliminate wrinkles and distortion during preparation of paraffin embedded tissue sections
- L.E.D. illuminated black background with transparent PYREX® dish; allows specimen to be easily seen
- Digital microprocessor control with advanced safety over temperature protection
- Supplied with:
- Removable PYREX® glass dish

Accessories

- PYREX® glass dish (part number: 904-0013)
- Certified to UL61010-1 and CSA C22.2



	Cardinal Health Deluxe Lighted Tissue Flotation Bath	Cardinal Health Lighted Tissue Flotation Bath
Material code	M7654-5A	M7654-1A
Overall dimensions (Width x Length x Height)	36.2 cm x 36.2 cm x 11.4 cm (14.25 in x 14.25 in x 4.5 in)	36.2 cm x 36.2 cm x 11.4 cm (14.25 in x 14.25 in x 4.5 in)
Dryer block	5.1 cm x 12.7 cm (2 in x 5 in)	N/A
Orienter block	5.1 cm x 5.1 cm (2 in x 2 in)	N/A
Net weight		5.45 kg (12 lb)
Shipping weight		6.8 kg (15 lb)
Accuracy		+/-1° C
Uniformity at 55° C		+/-1° C
Dryer block temperature	40 °C +/- 3 °C	N/A
Orienter block temperature	55 °C +/- 3 °C	N/A
Approvals	Certified to UL61010-01 and CSA C22.2	
Electrical	115 VAC 50/60 Hz 272 W	115 Vac 50/60 Hz 272 W

Cardinal Health Tissue Baths have a 2 year warranty on parts and workmanship.

Cardinal Health Standard Tissue Bath

Material code: M7655A

- Edge serves as drying surface for slides
- Surface is black porcelain enamel coated steel
- Equipped with on/ off switch, pilot light, thermostat and non-skid rubber feet
- Rapid heat-up

Supplied with:

- Bath cover
- Thermometer holder



Material code	M7655A
Overall dimensions (Diameter x Height)	27.9 cm x 9.2 cm (11 in x 3 ¾ in)
Chamber dimensions (Diameter x Height)	21 cm x 5.1 cm (8 ¼ in x 2 in)
Rim (width)	3.2 cm (1 ¼ in)
Net weight	2.3 kg (5 lb)
Shipping weight	2.7 kg (6 lb)
Capacity	1 L (1 inch from top)
Temperature range	35° C to 65° C
Electrical	115 VAC 50/60 Hz 190 W

Cardinal Health Paraffin Wax Dispenser

Material code: M7654-8

- Small compact footprint
- Heated valve ensures a consistent flow of molten wax
- Stainless Steel Outer housing
- Aluminum inner reservoir provided excellent thermal transfer and temperature control
- Dispenser valve allow precise control of flowing wax and can be locked in the open position if desired
- Inner reservoir Screen helps to prevent large particulate from clogging the dispenser valve
- Standard two year manufacturers warranty

Material code	M7654-8
Capacity	5 Liter (1.25 Gallon)
Temperature Range	40°C to 75°C +/- 2 °C
Approvals	Certified to UL 61010-1 and CSA C22.2
Electrical	115VAC 50/60 Hz 250 W



We are your lab partner

Expertise – Having the most experience meeting unique laboratory needs

Choice – Offering the most options in quality laboratory products, manufacturers and services

Efficiency – A comprehensive single resource capable of servicing your laboratory needs

Savings – Offering the lowest delivered cost for your laboratory products

Cardinal Health Slide Cabinet

Material code: M6360-2P

- Six cabinet drawers hold up to 5000 standard 1 inch x 3 inch microscope slides
- Cabinets can be interlocked up to 10 units high
- Plastic drawer dividers can be removed to store objects up to 5 cm (2 in) wide
- Drawers slide smoothly and easily
- Back of drawer acts as handle for transportation
- Durable cold rolled steel with baked on epoxy finish
- Supplied with sponge block spacers and Identification labels
- Additional spacers (ordered separately) Material code BO9010007B
- Cabinet Base with Rubber Feet (ordered separately) Material code M6362-1



For more information or to make a purchase contact your
Cardinal Health lab partner or

- Speak with a customer service representative at 800.964.5227
- Consult your 2008-2009 SP Lab Products Catalog
- Visit cardinalhealth.com



cardinalhealth.com

Cardinal Health
7000 Cardinal Place
Dublin, OH 43017