

Cardinal Health™ *H. pylori* Stool Antigen Test

Solution to detect active infection



Test. Treat. Retest to confirm eradication (no age restrictions)



Cardinal Health™ *H. pylori* Stool Rapid Test is a CLIA waived rapid immunoassay for the qualitative detection of *H. pylori* antigen from stool specimens collected in a clean air tight container, stored at 2 – 8° C until tested.

Ease of use

- Clear, objective results within 5 minutes
- Store specimen refrigerated but may be held up to 72 hours at 2 – 8° C prior to testing. If testing cannot be performed within this time frame, specimens should be frozen immediately on receipt and stored frozen ($\leq -20^{\circ}$ C) until tested


Precision

- Sensitivity: 91%
- Specificity: 92%

Value

- Accurate and rapid detection of active *H. pylori* infection
- No special equipment required

Serology testing is no longer recommended by AGA because it does not test for active infection and does not confirm for eradication.



Ordering information

Cat. No.	Description	Qty. Bx
B1077-25CH	Cardinal Health™ Rapid <i>H. pylori</i> Stool Antigen Test Kit	20 test/kit

Specimen Transport

CHB900312	Cardinal Health™ Clean or Empty Stool Transport	
CHB990125	Cardinal Health™ Stool Collection Device	

Accessories

CHV101-R	Cardinal Health™ Mini Vortexer, Red cup head	
CHV1000	Cardinal Health™ Vortex Mixer	

For more information or to order, contact your sales representative.