



Cardinal Health™ LIQUIBAND®
Topical Skin Adhesives

Trade your pen for
a new set of wings.

Traditional pen-style applicators have been the standard for years.

Now, there's a better way to work.


Our Cardinal Health™ LIQUIBAND® Topical Skin Adhesives feature a signature winged applicator to help improve safety for clinicians and patients.

For clinicians —


- Less chance of glove breach or sharps injury
 - With our applicator, the signature wings break the ampoule — not your fingers.
- One layer application in a single stroke

For patients —

- Provides water-resistant, microbial barrier
- Low exothermic reaction¹
 - Ideal for wound closure in the Emergency Department




Competing adhesives with accelerant




Low exothermic reaction

Cardinal Health™ LIQUIBAND®



Wound closure is an essential part of every OR procedure.

Some applicators clog so quickly that the adhesive flows for just a few minutes² and the rest dries in the applicator. Cardinal Health™ LIQUIBAND® Topical Skin Adhesives dry on the skin, not in the applicator.¹




Longer working time

- Applicator remains usable for up to 90 minutes after activation¹
- Adhesive flows easily so you can work at your own pace

With the competing pen-style applicator, the adhesive may not cover in a single stroke, requiring fine tuning such as filling in gaps, flattening ridges and eliminating bubbles.³ Many express the adhesive out all at once to avoid clogging and then brush a thick layer over the wound.

Only need **1** applicator for most procedures





Try the Cardinal Health™ LIQUIBAND® family of topical skin adhesives, including both Octyl and Butyl formulations in our signature winged applicator.



- Adhesive that polymerizes on skin, not in applicator¹
- Winged applicator that is safe, efficient, easy to activate and control¹

It simply works better — so you can too.

2-Octyl cyanoacrylate

	Cat No.	Description	Tip	Tint	Volume	Packaging
	CHLOT01-08	Cardinal Health™ Topical Skin Adhesive, LIQUIBAND® OCTYL	Dome	Violet	0.8mL	10 each/box, 6 box/case
	CHLOM01-04	Cardinal Health™ Topical Skin Adhesive, LIQUIBAND® OCTYL Mini	Dome	Violet	0.4mL	10 each/box, 6 box/case

n-Butyl cyanoacrylate

	Cat No.	Description	Tip	Tint	Volume	Packaging
	LQB003T	Cardinal Health™ Topical Skin Adhesive, LIQUIBAND® Butyl	Dome	Violet	0.5mL	12 each/box, 6 box/case
	LFC004T	Cardinal Health™ Topical Skin Adhesive, LIQUIBAND® Butyl Flow Control	Flow Control	Violet	0.5mL	12 each/box, 6 box/case

For more information, contact your local Cardinal Health sales representative or customer service at 800.964.5227.

¹ Data on file at Cardinal Health and Advanced Medical Solutions Ltd.

² Data on file at Cardinal Health and Advanced Medical Solutions Ltd. Ethicon Dermabond Advanced/Mini Directions for Use – Step 4

³ Key Group Research Study 2017 – data on file at Cardinal Health

⁴ Surgical Tissue Adhesive Wound Closure Length Study 2017 – data on file at Cardinal Health and Advanced Medical Solutions (Plymouth) Ltd. Compared to Ethicon wound coverage claim on webpage (15 cm). <http://www.ethicon.com/na/products/wound-closure/skin-adhesives/dermabond-advanced-topical-skin-adhesive>

LIQUIBAND is a registered trademark of Advanced Medical Solutions Ltd.

© 2019 Cardinal Health. All Rights Reserved. CARDINAL HEALTH, the Cardinal Health LOGO and ESSENTIAL TO CARE are trademarks of Cardinal Health and may be registered in the US and/or in other countries. All other marks are the property of their respective owners. Lit. No. 2MP19-905019 (5/2019)



Featuring



cardinalhealth.com/liquiband