

Cardinal Health™ Urinalysis System

Reduce variability. Reduce cost.



TRANSPORT RACK

Shown:
Catalog number
CH3048-2A

Urinalysis testing techniques can vary among different individuals. Reduce variability in routine urine examinations with the **Cardinal Health™ Urinalysis System**. Designed for consistent and cost-effective performance, our kits and components standardize critical steps in the macroscopic and microscopic examination of clinical urine specimens.

Cat. No.	Description	Qty
CH3050-1	Economy urine tubes 12mL, graduated, polystyrene	500 ea/cs
CH3050-2	Ultra urine tubes 12mL, graduated, break-resistant	500 ea/cs
CH3050-3	Caps for urine tubes	500 ea/cs
CH3050-4	Urine decanting pipettes	500 ea/pk, 10 pk/cs
CH3050-17	10-chambered slides	100 ea/pk, 10 pk/cs
CH3050-9	Ultra system pack for 1000 urinalysis determinations: 1000 ultra urine tubes 1000 urine decanting pipettes 100 10-chambered slides	1 kit/cs
CH3050-5	Ultra system pack with caps for 1000 urinalysis determinations: 1000 ultra urine tubes 1000 caps for urine tubes 1000 urine decanting pipettes 100 10-chambered slides	1 kit/cs
CH3048-2A	Complete system pack for urine collection/transport/testing up to 500 urine samples: 5 bags of 100 ultra urine tubes 5 bags of 100 urine cups 5 bags of 100 caps for urine tubes 5 bags of 100 labels 5 transport racks	1 kit/cs
CH2195-1	Complete system pack II for urine collection/transport/testing up to 500 urine samples: 5 bags of 100 economy urine tubes 5 bags of 100 urine cups 5 bags of 100 caps for urine tubes 5 bags of 100 labels 5 transport racks	1 kit/cs
CH3050-8	Economy system pack for 500 urinalysis determinations: 500 economy urine tubes 500 urine decanting pipettes 500 10-chambered slides	1 kit/cs

Designed for quality and convenience:

- Ready-to-use disposable components.
- Volume control for consistent microscopic examination.
- Components designed to work in conjunction with one another.
- Ultra urine tube (12mL) produced from special polystyrene to reduce the risk of breaking and cracking.
- Decanting pipette decants off the supernatant while maintaining 1mL of sediment at the bottom of the tube.
- Optically clear acrylic slide features an integrated coverslip with 10 chambers. Each chamber is individually sealed to prevent sample cross contamination.

We are your lab partner.

Expertise. Choice. Efficiency. Savings.

For more information or to order,

contact your Cardinal Health sales representative, or call **800.964.5227** or visit **cardinalhealth.com**.



cardinalhealth.com

Cardinal Health
Laboratory Products
Waukegan, Illinois