

The Unit Dose Pig provides radiation shielding. The top includes an o-ring to protect against possible fluid leaks.

The red disposal cap incorporates a biohazard symbol. The cap locks in place after the used syringe is returned to the insert, creating an OSHA-compliant point-of-use sharps container.

The single-use SECURE® Safety insert keeps the pig lead shield clean, helping to eliminate potential contamination of doses from previously used syringes.

The SECURE® Safety Insert is puncture resistant and leakproof on the top, bottom and sides.

The wider base and ridges minimize rolling.

