



COVIDIEN

positive results for life™

Uni-Patch™ **C Series (First Choice)** **Stimulating Electrodes**



C Series (First Choice)
Stimulating Electrodes

Long Lasting .040 inch Thick Hydrogel

- Provides reliable adhesion that conforms well over body contours.

Quick, Moderately Aggressive Tack

- Firm, yet gentle electrode to skin contact for multiple applications.

Highly Conductive Carbon Film

- Even current distribution.

Spun Lace Top Coat or Soft Foam Covers

- Flexible and conformable to a variety of application sites.

Soft Mold Flex Connector with Secure Lead Wire Adhesive Technology

- The connection socket easily adapts to variations in pin size and provides a secure connection.
- Wire adhesive is designed to anchor the lead wire into the pad and diminish separation.

Highly Conductive Pigtail Wire Strands

- Flexible and durable for multiple uses.

Double-Coated and Tinted Storage Liner

- User can return electrodes to either side of the storage liner.






Uni-Patch™ C Series (First Choice) Stimulating Electrodes

Ordering Information

First Choice Cloth Stimulating Electrodes

	Code	Model	Size
	EP85795	3100C	1.25" Round
	EP85765	3105C	2" Round
	EP85770	3110C	2.75" Round
	EP85775	3115C	2" Square
	EP85780	3120C	2" x 3.5" Rectangle
	EP85785	3125C	1.5" x 2.5" Oval
	EP85790	3130C	2" x 4" Oval

First Choice Foam Stimulating Electrodes

	Code	Model	Size
	EP85800	3150F	1.25" Round
	EP85805	3155F	2" Round
	EP85825	3160F	2.75" Round
	EP85835	3165F	2" Square
	EP85840	3170F	2" x 3.5" Rectangle

4 Electrodes per pack.

10-14 Applications.

Meets U.S., Canadian and European (CE) regulations. Made in USA.

www.covidien.com/uni-patch

COVIDIEN, COVIDIEN with logo, Covidien logo and *positive results for life* are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company. ©2015 Covidien.

H8823 – xM – 0215



1313 WEST GRANT BOULEVARD
WABASHA, MN
55981

800-328-9454
651-565-2601

WWW.COVIDIEN.COM