



**COVIDIEN**

*positive results for life™*

## Uni-Patch™ PolyHesive® Blue Gel Electrodes




**PolyHesive®  
Blue Gel Electrodes**



**The PolyHesive® Blue Gel  
Electrodes remove cleanly,  
painlessly and easily.**

Uni-Patch™ stimulating electrodes with proprietary PolyHesive® conductive blue gel are designed especially for patients with skin sensitivities. PolyHesive® bonds seamlessly to the skin and conforms to even minor irregularities.

PolyHesive® blue gel provides a clean, efficient, temporary bond to the surface of the skin. Extra thick gel promotes reliable adherence and reusability, while the high moisture content delivers excellent conductivity.



- Through air displacement, PolyHesive® sticks and conforms completely to skin irregularities, voids and pores.

# Uni-Patch™ PolyHesive® Blue Gel Electrodes

PolyHesive® hydrogel is a dependable choice for non-irritating, biocompatible skin adhesion.

## Ordering Information

	Code	Model	Size
	EP85205 (S-Series)	696SS	2" x 2" Square, 4/PK
	EP85215 (S-Series)	697SS	1.75" x 3.75" Rectangle, 4/PK
	EP85577 (ST-Series)	545-ST	2" x 2" Square, 4/PK
	EP85571 (ST-Series)	555-ST	1.75" x 3.75" Rectangle, 4/PK
	EP84604	639	4" x 6" Butterfly, 1/PK

[www.covidien.com/uni-patch](http://www.covidien.com/uni-patch)



COVIDIEN, COVIDIEN with logo, Covidien logo and *positive results for life* are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company. ©2015 Covidien.

H8962-1M-0715

1313 WEST GRANT BOULEVARD  
WABASHA, MN  
55981

800-328-9454  
651-565-2601

WWW.COVIDIEN.COM