



**COVIDIEN**

*positive results for life™*

## **Uni-Patch™ R Series (Re-Ply) Stimulating Electrodes**



**R Series (Re-Ply)  
Stimulating Electrodes**

### **Long Lasting .025 inch Thick Hydrogel**

- Provides reliable adhesion that conforms well over body contours.

### **Quick, Moderately Aggressive Tack**

- Firm, yet gentle electrode to skin contact for multiple applications.

### **Soft Mold Flex Connector with Secure Lead Wire Adhesive Technology**

- The connection socket easily adapts to variations in pin size and provides a secure connection.
- Wire adhesive is designed to anchor the lead wire into the pad and diminish separation.

### **Highly Conductive Pigtail Wire Strands**

- Flexible and durable for multiple uses.

### **Highly Conductive Carbon Film**

- Even current distribution.

### **Durable Tricot Covers**









- Helps to extend life of electrode.

### **Double-Coated and Tinted Storage Liner**

- User can return electrodes to either side of the storage liner.

# Uni-Patch™ R Series (Re-Ply) Stimulating Electrodes

## Ordering Information

	Code	Model	Size
	EP84770	650	1.5" x 2" Oval
	EP84820	652	2" Round
	EP84850	653	2.75" Round
	EP84910	654	2" x 2" Square
	ES82180	655	2" x 2" Square, Snap-Connection
	EP84772	657	2" x 4" Oval
	EP85040	658	2" x 4" Rectangle
	ES82200	659	2" x 4" Rectangle, Snap-Connection

4 Electrodes per pack.  
10-14 Applications.

Meets U.S., Canadian and European (CE) regulations. Made in USA.

[www.covidien.com/uni-patch](http://www.covidien.com/uni-patch)



COVIDIEN, COVIDIEN with logo, Covidien logo and *positive results for life* are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company. ©2015 Covidien.

H8824 – xM – 0215

1313 WEST GRANT BOULEVARD  
WABASHA, MN  
55981

800-328-9454  
651-565-2601

WWW.COVIDIEN.COM