



**COVIDIEN**

*positive results for life™*

## Uni-Patch™ S Series (Classic) Stimulating Electrodes

Classic electrodes feature our exclusive hydrogel with pure aloe vera. The thick hydrogel layer provides consistent adherence even after multiple uses.

### Silver Ink/Carbon Film Conductor

- Proven to efficiently disperse electrical current over the entire surface of the electrode.

### .040" Thick ALOE Hydrogel

- This gel is formulated with pure aloe vera, known for its hydrating properties and natural moisturizing effects.

### Intuitive Lift Tab

- Allows for easy removal of the electrode from the liner and from the patient.
- Helps to prevent accidental pulling on the lead wire, which can damage and shorten the useful life of the electrode. (Tab not available on TWF14 or TWF24)

### Soft Mold Flex Connector with Secure Lead Wire Adhesive Technology

- The connection socket easily adapts to variations in pin size and provides a secure connection.
- Wire adhesive is designed to anchor the lead wire into the pad and diminish separation.

### Soft Foam Backing

- Flexible and conformable foam provides comfort.

### Tinted Double Coated Storage Liner








- Visible liner makes electrode storage easy as either side of the liner can be used.



**S Series (Classic)  
Stimulating Electrodes**

# Uni-Patch™ S Series (Classic) Stimulating Electrodes

## Ordering Information

	Code	Model	Size
	EP84179	1504	0.75" x 1.375" Rectangle
	EP84159	2204	2" Square
	EP84169	2404	2" x 4" Rectangle
	EP84149	2604	2" Round
	EP85159	2804	2.75" Round
	EP84699	TWF14	1.5" x 2.5" Oval
	EP84599	TWF24	2" x 4" Oval

4 Electrodes per pack.  
10-14 Applications.

Meets U.S., Canadian and European (CE) regulations. Made in USA.

[www.covidien.com/uni-patch](http://www.covidien.com/uni-patch)

COVIDIEN, COVIDIEN with logo, Covidien logo and *positive results for life* are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company. ©2015 Covidien.

H8822 – xM – 0215



1313 WEST GRANT BOULEVARD  
WABASHA, MN  
55981

800-328-9454  
651-565-2601

WWW.COVIDIEN.COM